Poulsen & Graf

SUPERIOR LABORATORY PRODUCTS

OPTIFIX[®]

Bedienungsanleitung
Operating manual
Mode d'emploi
Istruzioni per l'uso
Instrucciones de
manejo
操作说明

FORTUNA[®]

Table of Chemical Resistance

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Preface

We thank you for the purchase of an OPTIFIX® bottle top dispenser. You made a good choice. With 40 years of experience in making dispensers we guarantee a high standard of product quality and safety. In order to guarantee the quality and safety of the dispenser for long term dosing of liquids, correct operation and maintenance are required.

Therefore please read this operating manual carefully before using the dispenser and take note of the warning instructions as well as the corresponding regulations. Please keep it for later reference.

The pack contains: 1 dispenser (valve block unit with finger guard), 1 discharge tube set, 1 intake tube, 2 adapters GL 40+45 (model **BASIC** and **SOLVENT**: made of PP; model **SAFETY**, **SAFETY S** and **HF**: made of PTFE), 1 operating manual. Please check this before using the dispenser.

Please pay attention to the colour code of our bottle top dispensers:

■ **BASIC**, ■ **SOLVENT**, ■ **SAFETY**, **SAFETY S**, ■ **HF**

Operating Manual Bottle Top Dispenser OPTIFIX®

I. Safety Instructions

1. Before using the dispenser please read the operating manual carefully and follow the safety instructions. ⚠ **This operating manual cannot reasonably mention all safety risks that might be associated with the use of the dispenser. It is therefore in the responsibility of the user to consult where necessary and establish appropriate health and safety practices to determine the applicability of regulatory limitations prior to use.**
2. Observe general notes on hazards and safety regulations, e. g. wear protective clothing, gloves and glasses.
3. Observe any hazard symbols on reagent bottles.
4. Use only appropriate containers to dispense liquids. If using small bottles, please use a stand.
5. Always make sure that the discharge tube never points at the user or any other person.
6. When transporting the dispenser please ensure that the piston is in zero-position and no liquid is in the cylinder.

The stopper cap must close the discharge tube and the volume setting system must be set at zero-position! Never carry the dispenser by the finger guard, only by the PTFE-valve block!

7. Before use, please check correct assembly of the dispenser with the bottle, that the bottle is properly fitted and closed to full lock. This will be possible using standard glass laboratory bottles with a standard thread (e.g. Schott bottles). By use of short neck bottles, bottles or containers with a short thread, or other bottles where full tightening of the bottle top dispenser onto the bottle cannot be achieved, the dispenser should be screwed onto the bottle using an adaptor to ensure secure operation. Most adaptors are readily available from our catalogue or can be supplied according to special requirements on request. **The stopper cap of the discharge tube must never obstruct the dosing jet, in order to avoid dangerous splashing.**

8. To guarantee trouble-free working please clean and use the dispenser according to the operating manual.

9. Repairs if necessary must be carried out only by the manufacturer or authorized service stations using original parts.

10.  For safety reasons we recommend use of adaptors made of PTFE which are available as an option. The models **SAFETY**, **SAFETY S** and **HF** must be used with adaptors made of PTFE.

11. To minimize contamination from the sapphire valve balls before dispensing of hydrofluoric acid with the HF, we recommend discarding approx. 2 ml.

II. The World of OPTIFIX®

OPTIFIX® BASIC* (Art. No. 101 080)

Bottle top dispenser for aqueous liquids, acids, bases and crystallizing liquids. Piston made of glass, solid PTFE coated, patented. Safety dispensing cylinder, plastic coated (except 200-500 ml). Valve block made of pure PTFE. Valve system without metal springs. FIX adaptor.

OPTIFIX® SOLVENT* (Art. No. 101 081)

Bottle top dispenser for solvents. Piston made of precision glass. Safety dispensing cylinder, plastic coated (except 200-500 ml). Valve block made of pure PTFE. Valve system without metal springs. FIX adaptor.

FORTUNA UNIVERSAL® OPTIFIX® HF (Art. No. 101 106)

Bottle top dispenser for extremely hazardous chemicals, especially hydrofluoric acid. Piston made of glass, solid PTFE coated (patented), with sealing lip. Dispensing cylinder and valve block made of pure PTFE. Valve system without metal springs. Ventilation with ball valve. FIX adaptor.

OPTIFIX® SAFETY* (Art. No. 101 107)

Bottle top dispenser for strong acids or hazardous chemicals. Piston made of glass, solid PTFE coated (patented), with sealing lip. Safety dispensing cylinder, plastic coated. Valve block made of pure PTFE. Valve system without metal springs. FIX adaptor made of solid PTFE.

OPTIFIX® SAFETY S* (Art. No. 101 107S)

Bottle top dispenser for strong or fuming acids or highly odorous chemicals. Piston made of glass, solid PTFE coated (patented), with sealing lip. Safety dispensing cylinder, plastic coated. Valve block made of pure PTFE. Valve system without metal springs. Ventilation with ball valve. FIX adaptor made of solid PTFE.

 * Warning: not suitable for hydrofluoric acid

Never use the OPTIFIX® Dispensers with:

- liquids attacking Al_2O_3 and PTFE
- liquids attacking borosilicate glass (except with the **HF**)

Warning:

The **BASIC** and **SOLVENT** model come with adaptors made of polypropylene to fit on most popular containers. However, when dosing extremely hazardous chemicals, for safety reasons it is recommended to use adaptors made of PTFE for **BASIC** and **SOLVENT** model which are available as an option.

III. Modern Technique

Finger guard made of PP. All descriptions on the finger guard are printed and the graduation on the scale is lasered in order to guarantee permanence.

Volume setting system Easy-to-use volume setting system with fine adjustment for calibration. Suitable for use by left- and right-handers.

Dispensing cylinder Glass cylinder, vacuum-calibrated, plastic coated (200-500 ml without coating), with stopper ring. Black graduation. **HF** with dispensing cylinder made of pure PTFE.

Valve block made of pure solid PTFE, with FIX Adaptor.

Valve seats made of glass, precision ground for precision sapphire valve balls. No metal springs. Works by means of gravity. **HF** with valve seats made of sapphire.

Dispensing piston Model **BASIC**, **SAFETY**, **SAFETY S** and **HF** glass piston, coated with solid PTFE, patented. Sealing lip with model **SAFETY**, **SAFETY S** and **HF**. Model **SOLVENT** with glass piston.

Discharge tube set consists of: discharge tube, support, screw cap, thread adaptor and stopper cap.

IV. Handling


1. Assembling / Filling

Before using the dispenser please check the dispenser for any damage or glass breakage.

Push intake tube over the hose connection of the PTFE valve block.

Screw dispenser direct onto containers with thread GL 32 or use adaptors to fix on any other suitable threads or ground joint necks. Check the correct fitting. The dispenser can be turned on the flask and fixed in any direction by using the built-in FIX adaptor.

Screw the discharge tube set onto the threaded nipple of the valve block by hand. **Caution:** do not use any tools, otherwise the valve system may be damaged. Remove stopper cap from the discharge tube.

⚠ Warning:

If the stopper cap is not removed when dispensing, there may be an overpressure in the channel system and uncontrolled dispensing may be the result.


To fill the dispenser the first time lift and lower piston a short distance several times until no air bubbles are seen in the cylinder.

HF: Lift and dose the piston several times to fill the channel system.

2. Volume setting

Loosen the nut of the volume setting system, set the desired volume. The upper edge of the fine adjustment screw shows the volume set. Fix the screw.

For calibrating purposes turn the screw of the fine adjustment. Turn downwards to increase the volume, turn upwards to decrease the volume.


3. Dispensing procedure

a) For safety make sure that the discharge tube does not point at the user or any other person. By using the built-in FIX adapter it is possible to fit the dispenser as required.

Remove stopper cap by holding the discharge tube tight with other hand.

b) Set the desired volume with the volume setting system.

c) Lift the dispensing piston up with the finger guard. Do not touch the piston itself to avoid any unlikely contact with chemicals. The piston stops at the volume set. The graduation on the glass cylinder is an additional optical control of the volume set or air-bubbles. There is no graduation on the PTFE-cylinder of the **HF**.

d) The models **BASIC** and **SOLVENT** work by gravity due to its design. To dispense liquid let the piston go; it will fall by its own weight. If this is not the case, please apply the dosing piston manually.

The models **SAFETY**, **SAFETY S** and **HF** are equipped with a special piston seal. To dispense liquid push the piston down to the zero-position using the finger guard.

⚠ The stopper cap of the discharge tube may never come into the dosing jet, in order to avoid dangerous squirting.

e) Avoid uneven and strong piston movements or jerking of the stopper ring.

f) Replace stopper cap after dispensing by holding the discharge tube tight with other hand.

Warning:

The piston of the models **BASIC** and **SOLVENT** run down to the zero-position automatically under their own weight due to its design. If this is not the case, please apply the dosing piston manually. With the models **SAFETY**, **SAFETY S** and **HF** with special seal the piston must be pushed down manually. Therefore chemicals may remain in the dispensing cylinder. To avoid unintentional dispensing, stopper cap must be placed on the discharge tube after every use. When dispensing observe safety precautions!

4. Change the dosing medium bottle

To change the dosing medium bottle, please proceed as follows.

1. move the dosing piston manually all the way down and make sure that there is no more liquid in the dosing cylinder.
2. loosen the fixed adapter of the dispenser and lift the dispenser together with the tube until the tube is no longer in the liquid.
3. dispense the remaining liquid in the dispenser by lifting and lowering the dispensing piston once or twice. Then remove the suction hose.

 **Warning!** Wear suitable protective clothing. Residual dispensing liquid may drip from the intake tube. Take appropriate precautions. When changing the dosing liquid, it is recommended to clean the unit.

V. Operating Temperature

The operating temperature of the OPTIFIX Dispensers is from +15° C to +35° C. (storage temperature from -20° C to +50° C). Under consideration of the specified precision you have to take care that the operating temperature is between +20° and +30° C.

VI. Sterilization

Model	autoclavable +121°C (248°F) 1 bar	not autoclavable	autoclavable on request
BASIC (2-100 ml) ⁽¹⁾	X	-	-
SOLVENT (2-100 ml) ⁽²⁾	X	-	-
BASIC (200, 300, 500 ml)	-	-	X
SOLVENT (200,300,500 ml)	X	-	-
SAFETY/SAFETY S/HF ⁽³⁾	-	X	-

- ⁽¹⁾ The **BASIC** piston has to be sterilized separately. Before sterilisation loosen screw of the volume setting system. **The dispensing piston has to be cooled down to room temperature (after approx. 2 hours) before reassembling.**
- ⁽²⁾ The dispensing piston of the **SOLVENT** can remain in the cylinder during the sterilization.
- ⁽³⁾ Please clean with appropriate cleaning solutions or ultrasound. Sterilisation of **SAFETY**, **SAFETY S** and **HF** OPTIFIX® dispensers will put excessive stress on the dispensers both thermally and mechanically. As a consequence sterilisation of OPTIFIX dispensers will invalidate the warranty offered on the dispensers.

VII. Transport and Storage

1.  When moving the dispenser never carry by the finger guard but always by the PTFE- valve block to avoid unintentional dispensing of liquid! Replace the stopper cap before transporting and fix the volume setting system at the zero-position.
2. Please transport and store the dispenser only in the original box to avoid damage. Storage temperature: -10° C (+14° F) to +50° C (+122° F).

VIII. Maintenance and Cleaning

1. To guarantee trouble-free use for many years it is necessary to clean the dispenser regularly (e.g. weekly before a longer dosing break). If there is any sign of a potential malfunction (e.g. piston is difficult to move) never use force. Immediate cleaning should be carried out.

2. Do not leave the dispenser empty without cleaning to avoid clogging the valve balls! If this does happen soak the dispenser for some hours in cold water or clean it using ultrasound. Once the piston can be moved, dismantle and remove any residual impurities by cleaning with a normal solution.

Disassembly of the Dispenser:

a) Remove the intake tube.

b) Unscrew the discharge tube (caution- do not remove the thread nipple from the valve block)

c) Pull the finger protection cap up to the top, loosen the volume setting and remove it by pulling it down

d) Pull the finger protection cap and dispensing piston out of the dispensing piston

3. When changing the dispensing liquid the dispenser should be cleaned.

Please note: When mounting the SOLVENT Dispenser the glass piston has to be screwed into the glass cylinder carefully otherwise the piston and cylinder can get jammed.

IX. Repair Service

1. All OPTIFIX® dispensers are made from highest quality materials and are almost maintenance-free and reliable. If a repair should be necessary, it must be done by the manufacturer or an authorized service station. This does not refer to parts available as spare parts (finger guard with piston, volume setting system, discharge tube set, tubes, adaptors).

2. Before returning the dispenser please do not forget to register it with the manufacturer in advance. Make sure that it is perfectly cleaned to avoid any contamination. Please also enclose a filled in safety declaration.

3. Please use the dispenser following to the operating manual to be entitled to the warranty. Repairs must be carried out by the manufacturer or authorized service stations.

4. Please note that the HF dispensing cylinders and dispensing pistons are not interchangeable.

X. Disposal

When disposing the dispensers please pay attention to the national disposal regulations.

XI. Technical Data (EN ISO 8655/5)

1. Volume Range - Subdivisions - Accuracy

Volume range in ml	Graduation in ml	Accuracy A in %	Precision CV in %
0,5 - 2,0	0,1	$\leq \pm 0,6$	$\leq \pm 0,2$
1,0 - 5,0	0,1	$\leq \pm 0,6$	$\leq \pm 0,2$
2,0 - 10,0	0,2	$\leq \pm 0,6$	$\leq \pm 0,2$
5,0 - 30,0	0,5	$\leq \pm 0,6$	$\leq \pm 0,2$
10,0 - 50,0	1,0	$\leq \pm 0,6$	$\leq \pm 0,2$
20,0 - 100,0	2,0	$\leq \pm 0,6$	$\leq \pm 0,2$
40,0 - 200,0	5,0	$\leq \pm 0,6$	$\leq \pm 0,2$
60,0 - 300,0 *	5,0	$\leq \pm 0,6$	$\leq \pm 0,2$
100,0 - 500,0 *	5,0	$\leq \pm 0,6$	$\leq \pm 0,2$

*not acc. to ISO 8655/5

The OPTIFIX® dispensers are calibrated to deliver (EX).

The data for accuracy and precision are established according to the German calibration ordinance (PTB) at a room temperature of +20° C (+68° F) by using distilled water. When using chemicals with a different viscosity or density recalibration on an analytical balance is needed.

2. Packing

OPTIFIX® dispenser 2ml, 5ml, 10ml, 30 ml, 50ml and 100ml:
single carton, 29 x 16,5 x 8 cm.

OPTIFIX® dispenser 200ml, 300ml and 500ml: single carton,
46,2 x 18,5 x 14,5 cm.

3. Weight, Dimensions

Volume range in ml	Weight in g appr.	Height in mm
0,5 - 2,0	225	150
1,0 - 5,0	240	165
2,0 - 10,0	300	180
5,0 - 30,0	485	235
10,0 - 50,0	605	245
20,0 - 100,0	800	275
40,0 - 200,0	1500	365
60,0 - 300,0	1800	380
100,0 - 500,0	1900	435

XII. Equipment

All models come in a break-proof single carton with 1 intake tube and 1 discharge tube set. The models **BASIC** and **SOLVENT** are additionally equipped with 2 PP adapters, the models **SAFETY**, **SAFETY S** and **HF** are equipped with 2 PTFE adapters.

XIII. Warranty

The OPTIFIX® dispenser is probably the most reliable bottle top dispenser in the world. The manufacturer guarantees manufacture and original parts for a 12 month period but excluding all claims for damages, conversion and other claims. Under this warranty all parts that are defective due to faulty materials or production will be replaced. The manufacturer is not responsible for any damages or problems caused by transport or improper handling. All claims must be reported to the supplier. After a warranty repair the warranty period will not be extended. This warranty is not transferrable and void if repairs are carried out by unauthorized service stations.

Description of the malfunction	Possible reason	Solution
Air bubbles in the dispensing cylinder	The dispensing piston has been pulled up too quickly. The system has not been vented.	Slowly and uniformly pulling out and pressing down the dispensing piston. Repeated filling and emptying of the dispensing cylinder
The dispensed volume is too low.	The dispensing piston has been pulled up too quickly. The system has not been vented. There are impurities in the intake tube.	Slowly and uniformly pulling out and pressing down the dispensing piston. Repeated filling and emptying of the dispensing cylinder Remove and clean the intake tube
The dispensing piston is sluggish.	There is contamination between the dispensing piston and the dispensing cylinder. The PTFE piston is swollen.	Disassembly and cleaning of dispensing piston and dispensing cylinder Disassembly of the dispenser, put the PTFE piston together with the finger protection cap in the freezer for a short time. Check the dosing liquid, liquids with Solvent-characteristics can only be dispensed with our SOLVENT
Fluid is leaking between the valve block and the discharge tube	The discharge tube has not been installed correctly.	Unscrew the discharge tube and screw it again. Do not overtighten!
The dispensing cylinder can not be filled.	The volume setting system sits at the top. The intake tube is missing.	Adjust the volume setting system to the desired volume Unscrew the dispenser from the bottle and insert the intake tube
Dispensing is not possible.	The discharge tube is blocked.	Unscrew the discharge tube and clean it.

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All rights reserved, especially the rights to illustrations and original texts in whole or in part, for any and all reproduction.
OPTIFIX Manual Rev. 10/2025, Copyright by Poulsen & Graf GmbH,