General Catalogue

Poulten & Graf

SUPERIOR LABORATORY PRODUCTS

Catalogue Edition 2022-06

Copyright, Layout, Text, Photos

Poulten & Graf GmbH,

Photo Studio Eilers, Stadtprozelten

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Copyright, Layout, Text, Fotos

Poulten & Graf GmbH,

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Experience perfection
through a combination
of quality, precision
and design!

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02. Liquid Handling · Dosiertechnik

No matter whether it is manual or automatic, aqueous solution or aggressive acid - we have the right device for your application!

Due to our long experience in the production of manual bottle top dispensers and automatic dosing systems, our products are easy to use, safe and precise. They are indispensable in everyday life in the laboratory.

We are proud to present our wide product range, which extends from bottle top dispensers and microliter pipettors to an automatic dosing system. These cover most laboratory applications and a number of chemical and pharmaceutical manufacturing uses as well. With one of the widest product range on the laboratory market we offer real choice and expertise to our customers.

We can answer your dosing questions. Please contact us!

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	Acids Säuren		HF HF	Alkaline Sol. Laugen	Aqueous Sol. wässrige Lösungen	Solvents Lösungsmittel
	weak	strong				
	schwach	stark				
 FORTUNA® POLYFIX® 101 070*	not recommended	not recommended	not recommended	recommended	recommended	preferred recommended
 FORTUNA® POLYFIX® 101 071*	recommended	recommended	not recommended	preferred recommended	preferred recommended	recommended
 FORTUNA® POLYFIX® 101 072*	recommended	recommended	not recommended	recommended	recommended	preferred recommended
 FORTUNA® OPTIFIX® BASIC	preferred recommended	recommended	not recommended	preferred recommended	preferred recommended	recommended
 FORTUNA® OPTIFIX® SOLVENT	recommended	recommended	not recommended	recommended	recommended	preferred recommended
 FORTUNA® UNIVERSAL® OPTIFIX® HF	preferred recommended	preferred recommended	preferred recommended	preferred recommended	preferred recommended	recommended
 FORTUNA® OPTIFIX® SAFETY	preferred recommended	preferred recommended	not recommended	preferred recommended	preferred recommended	recommended
 FORTUNA® OPTIFIX® SAFETY S	preferred recommended	preferred recommended	not recommended	preferred recommended	preferred recommended	recommended
 OPTIMAT® Pumps Pumpen	preferred recommended	preferred recommended	not recommended	preferred recommended	preferred recommended	recommended

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* all media against which polypropylene is resistant

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For all other liquids please contact us!

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FORTUNA® POLYFIX® Dispenser

- FORTUNA® POLYFIX® - with glass piston**
- FORTUNA® POLYFIX® - with PTFE-coated piston**
- FORTUNA® POLYFIX® - with amber glass cylinder and glass piston**

- visible dispensing procedure
- easy-to-use volume setting system
- valve block made of polypropylene with FIX adapter
- according to EN ISO 8655-5

Due to the fact that PP is suitable for a lot of chemicals, the POLYFIX® dispenser can be used for most liquids.

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according to EN ISO 8655-5

101 070

- with glass piston
- suitable for all chemicals (incl. solvents), that do not attack polypropylene and don't crystallize

101 071

- with coated PTFE-piston
- suitable for acids and alkaline solutions that do not attack polypropylene

101 072

- with amber glass cylinder and glass piston
- suitable for all light sensitive liquids that do not attack polypropylene

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comes with 1 PTFE intake tube, 2 PP adaptors (GL 40 + 45), packing unit: 1

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Replacement Parts - for 101 070

101 077

Glass piston, complete with finger guard and volume setting system

101 075

Valve block unit, complete with dispensing cylinder and discharge tube

- for 101 071

101 078

PTFE coated piston, complete with finger guard and volume setting system

101 075

Valve block unit, complete with dispensing cylinder and discharge tube

- for 101 072

101 077

Glass piston, complete with finger guard and volume setting system

101 079

Valve block unit, complete with amber dispensing cylinder and discharge tube

- for 101 070 - 101 072

101 089

Stopper caps for discharge tubes, clear

101 09338

Intake tubes made of PTFE (350 mm; 6/7 mm diameter), pack 2

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Replacement Parts for POLYFIX® Dispenser

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FORTUNA® OPTIFIX® Dispenser

Advantages at a glance

- Best safety by use of high grade material
- Visible dosing procedure
- No contamination by metal parts, e.g. springs
- High accuracy, according to EN ISO 8655-5
- Premium quality
- BASIC and SOLVENT - wide volume range up to 500 ml

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02. Liquid Handling · Dosiertechnik

Dispensing Piston

- proven technology
- made of glass, solid PTFE-encased, exception: SOLVENT
- virtually indestructible

Precision Dispensing Cylinder

- made of glass, dosing procedure visible (except type HF, made of PTFE)
- coated with plastic - protection against breakage (expect 200 - 500 ml)
 - graduation clear readable
- easy to dismantle and to clean

Valve Block

- 100 % tightness
- no metal springs - no contamination of the dosing liquid

Finger Protection Cap

- ergonomic type
- timeless design
- with colour-code-system

Volume Setting System

- fast and easy volume setting
- additional fine adjustment
 - for right- and left hand use

FIX Adaptor

- to align the dispenser on a bottle
- made of PP for BASIC + SOLVENT
- made of PTFE for SAFETY +HF

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according to EN ISO 8655-5

101 080

**for dispensing light acids,
bases, aqueous solutions
and crystallizing liquids**

- with solid PTFE-encased dispensing piston
- dispensing cylinder made of glass
- visible dosing procedure
- no metal spring, no contamination of the dosing liquid
- volume up to 500 ml

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The bottle is not included in delivery (Matching bottles you find on p. 182 ff)

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101 081

for dispensing solvents

- with glass piston
- dispensing cylinder made of glass
- visible dosing procedure
- no metal spring, no contamination of the dosing liquid
- volume up to 500 ml

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101 081
200 - 500 ml

comes with 1 PTFE
intake tube, 2 PP
adaptors (GL 40 + 45),
packing unit: 1

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according to EN ISO 8655-5

101 106

For all liquids that can be handled with the OPTIFIX dispensers and for hydrofluoric acid

- made of PTFE
- with additional safety valve in the ventilation tube
- with PTFE-encased dispensing piston and PTFE-cylinder
- FIX adaptor made of solid PTFE
- no metal spring, no contamination of the dosing liquid
- volume up to 50 ml

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comes with 1 PTFE intake tube, 2 PTFE adaptors (GL 40 + 45), packing unit: 1

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according to EN ISO 8655-5

101 107

for dispensing strong acids and aggressive chemicals

- PTFE-encased dispensing piston with sealing lip
- dispensing cylinder made of glass
- visible dosing procedure
- FIX adaptor made of solid PTFE
- no metal spring, no contamination of the dosing liquid
- volume up to 100 ml

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according to EN ISO 8655-5

101 107 S

for dispensing strong fuming acids and very abrasive liquids

- with safety stopcock and additional safety valve in the ventilation tube
- PTFE-encased dispensing piston with sealing lip
- dispensing cylinder made of glass
- visible dosing procedure
- FIX adaptor made of solid PTFE
- no metal spring, no contamination of the dosing liquid
- volume up to 100 ml

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Replacement Parts - for 101 080

101 000

Piston BASIC, solid PTFE encased, with finger guard and volume setting system

101 004

Valve block unit with dispensing cylinder

- for 101 081

101 001

Piston SOLVENT, made of glass, with finger guard and volume setting system

101 004

Valve block unit with dispensing cylinder

- for 101 107

101 003

Piston SAFETY, solid PTFE encased, with sealing lip, finger guard and volume setting system

101 004

Valve block unit with dispensing cylinder

- for 101 107S

101 003

Piston SAFETY S, solid PTFE encased, with sealing lip, finger guard and volume setting system

101 004S

Valve block unit with dispensing cylinder, with safety valve

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Replacement Parts for OPTIFIX® Dispenser

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Volume setting system

101 000KS

Volume setting system for BASIC, SOLVENT, HF, SAFETY, SAFETY S, POLYFIX, 2-100 ml

101 000KM

Volume setting system for BASIC, SOLVENT, 200 - 500 ml

Intake tubes

101 09336

Intake tubes for HF, 5 ml, made of PTFE

101 09338/48/51/58

Intake tubes for BASIC, SOLVENT, HF (10+30 ml), SAFETY, SAFETY S, POLYFIX, made of PTFE

101433

Tubing, made of PTFE, sold in varying lengths (state length with order)

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Discharge tube sets

101 089

Stopper caps for discharge tubes, natural (also for POLYFIX®)

101 094

Discharge tube set (tube with stabiliser, fixing screw, cap)

101 095

Discharge tube set with Safety-Stopcock (tube with stabiliser, 2 fixing screws, PTFE-stopcock with dispensing tip), completely interchangeable; prevents dripping after the dispensing procedure

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made of PTFE

101 105

Humidity absorber, made of PTFE, for hygroscopic liquids, single pack

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101 104

for OPTIFIX® SAFETY,
SAFETY S, OPTIMAT® 3

- flexible remote dispensing
- suitable for easy serial dispensing
- tubing and dispensing valve made of PTFE
- with new metal spring, chemical resistant
- prevents liquid from coming out and stops liquid flowing back

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**The WANDERER
for flexible
remote dispensing**

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Packing unit: 1

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FORTUNA® OPTIMAT® 3

Automatic programmable Dispensing Station

Technical Features

- Precision drive
- Colour touch screen
- Suitable for continuous operation
- 7 pump sizes available up to 300 ml

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New Features

- Large touch screen display for better reading precision
- New software for better handling
- on request language selection possible

Benefits

- designed to give years of continuous day-in day-out service
- gives precision at the highest stage
- dispenses even difficult and aggressive liquids
- extremely robust
- dispensing without corrosion or contamination

Typical Applications

- calibration of volumetric apparatus
- extremely accurate dispensing of strong acids, alkalis, high value chemicals in the pharmaceutical and cosmetic industries

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FORTUNA® OPTIMAT® 3

Automatic Dispensing Station

**suitable for 1 dispensing
pump, with RS 232
connection interface**

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Art.-No.
Art.-Nr.

150 0003

**Pump and tubes are
not included.
Please order separately!**

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FORTUNA® OPTIMAT® 3 dispensing pump

101 126

for dispensing strong fuming acids and very abrasive liquids, single pack (intake and discharge tube sets are not included)

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Pumps for OPTIMAT® 3 dispensing station:

Packing unit: 1

Pumpen für OPTIMAT® 3 Dosierstation:

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
101 126

Volume Volumen	Grad. Grad.	A* _{≤ ±} R* _{≤ ±}	CV* _{≤ ±} V* _{≤ ±}	max.V max.V	Art.-No. Art.-Nr.
0 - 2 ml	0.1 ml	± 0.6 %	± 0.2 %	1,90 ml/s	101 12627
0 - 10 ml	0.2 ml	± 0.6 %	± 0.2 %	3,35 ml/s	101 12638
0 - 30 ml	0.5 ml	± 0.6 %	± 0.2 %	9,60 ml/s	101 12644
0 - 50 ml	1.0 ml	± 0.6 %	± 0.2 %	15,20 ml/s	101 12648
0-100 ml	2.0 ml	± 0.6 %	± 0.2 %	24,30 ml/s	101 12651
0-200 ml ⁽¹⁾	5.0 ml	± 0.6 %	± 0.2 %	31,30 ml/s	101 12653
0-300 ml ⁽¹⁾	5.0 ml	± 0.6 %	± 0.2 %	31,80 ml/s	101 12656

* A=Accuracy, CV=Coefficient of variation *R=Richtigkeit, V= Variationskoeffizient
(¹) additional to ISO 8655 · in Ergänzung zur ISO 8655

101 132

101 134

Accessories for OPTIMAT® 3 and OPTIMAT® 3 dispensing pumps

101 132
Intake tube
length 1 m

101 134
Discharge tube set
made of PTFE, length 0.75
or 1.50 m

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155 00040
Serielles Verbindungskabel

**Tubes for
150 0003 OPTIMAT® 3
dispensing station:**

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Packing unit: 1

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110 325

With the desalination unit LAB-IoN L 2 demineralised water can be produced in the simplest way.

Advantages at a glance:

- output approx. 30 - 40 l/h
- used cartridges are simply replaced
- connection tube to a 3/4" tap

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Desalination unit LAB-IoN L 2 comes with:

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Packing unit: 1

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106 90003

The Water Distillator produces distilled water from tap water. The mono distillate produced is free of germs and dissolved salts. It complies with DAB 8 regulations.

106 900

Compact design, an automatic level regulation system, easy cleaning and a temperature sensor that switches off the distilling system in case of trouble, all argue for the Water Distillator.

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Precise and reliable –
we know what is
important to you.

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Poulten & Graf Volumetric Glassware - tested precision-made glassware

Poulten & Graf offers a wide range of Volumetric Glassware under the trade mark  FORTUNA®.

The whole range of Volumetric Glassware made by Poulten & Graf is produced according to British BS, European EN or international ISO standards.

There is also a specialised range made according to US Pharmacopoeia (USP) standards.

Poulten & Graf Volumetric Glassware is produced according to the highest quality standards. Control of raw materials, production and quality work hand in hand and guarantee the high quality of the final product.

Poulten & Graf Volumetric Glassware, conformity certified, class A, is produced with a LOT number and can be delivered with a batch certificate (except Original FORTUNA® graduated and bulb pipettes with suction piston are delivered without certificate).

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Graduation

We distinguish between two different graduation types:

- for Class A or AS products: complete ring marks
- for Class B products: short graduations

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Calibration

Generally we differentiate between volumetric glassware which is calibrated to contain (IN) or calibrated to deliver (EX). Flasks, measuring and mixing cylinders are calibrated to contain (IN) so the liquid poured out corresponds with the volume indicated. Pipettes and burettes, however, are calibrated to deliver (EX). Provided the required waiting time is observed, the delivered liquid corresponds with the volume indication.

 FORTUNA® glassware is individually calibrated (no mass production technique and no approximations).

Special Note: When users wish to check the accuracy of volumetric glassware, care should be taken not only to avoid parallax error, but also to avoid a false view of the meniscus. Use of a black ring or a clip approx. 1 mm below the mark will give a black reflection on the bottom of the meniscus and will dramatically improve the accuracy of the reading (see ISO 4787).

We can supply suitable black clips on request.

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Colour marking

The colour markings on **volac** FORTUNA® volumetric glassware are applied by silk screen printing. We only use diffusion or enamel colours. The diffusion colour - using amber stain - is completely absorbed into the glass surface during the firing process. It cannot be removed except by processes which remove the glass surface itself. It is therefore resistant to cleansing and even to very aggressive chemicals.

All other colours - e.g. blue, white, yellow, green - are enamel or melting colours, which are melted onto the glass surface at high temperature.

volac FORTUNA® glassware is individually printed. Particularly for graduated glass (graduated pipettes, burettes, cylinders) individual printing is the only way to ensure that the top and bottom of the scale match individual calibration exactly.

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Classes

In accordance with the ISO standard: FORTUNA®

Class	Certificate	Individually Numbered	Products
B	no	no	cylinders
A/AS	download batch certificate	on request	cylinders, flasks, burettes, graduated and bulb pipettes

In accordance with the ISO and USP standards: FORTUNA®

Class	Certificate	Individually Numbered	Products
A/AS	download batch certificate	no	burettes, graduated and bulb pipettes, cylinders, flasks
A/AS + WAC	single certificate	yes	cylinders, flasks
USP	single certificate	yes	burettes, graduated and bulb pipettes, cylinders, flasks

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Why use volumetric glassware according to USP-standard?

In the United States and in Arabic Countries, as well as in many Asian Countries there is strong demand for volumetric glassware according to USP- standards. Essentially, the tolerances of measuring devices made to these standards differ to a great extent from the European norm, ISO, and are much tighter.

If a company produces goods that have to comply with the American FDA import controls or come within their responsibility, they have to meet certain regulations (GMP).

GMP stands for Good Manufacturing Practice. This is a set of regulations, which ensures good quality, good production processes and appropriate production environments in the production of pharmaceutical products and chemical agents, as well as cosmetics, food and animal feed.

Amongst other things the focus lies in a complete and closed traceability of the products. In addition, the traceability and control of the manufacturing process (tracking), plays an important role.

To ensure the tracking of the volumetric glassware which Poulten & Graf produces, individual serial numbers and quality certificates are provided. Manufacturers, who do not follow the GMP and cannot ensure complete quality control, face warnings and penalties. As a last resort, producers, who aim to sell their products mainly in territory controlled by the FDA, can lose their permission to sell into the US market. To ensure maximum reliability, precision and traceability, choose **volac** FORTUNA® volumetric glassware, produced according to USP- standard, with individual serial numbers and certificates.

volac FORTUNA® Volumetric Glassware Specifications

Volumetric Flasks Class A	USP is tighter than ISO
Bulb Pipettes Class A	some USP specifications are the same as ISO but some are tighter
Graduated Pipettes Class A	again some USP specifications are tighter than ISO but some are the same
Measuring Cylinders Class A	USP specifications are unique
Burettes Class A	USP and ISO are identical

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FORTUNA® Burettes

 FORTUNA® burettes and automatic burettes are produced to the highest quality standards and only according to **Class AS (EN ISO 385)**.

Reading: With ring marks and Schellbach stripes, precise working is possible, reading the volume at the touch point of the two tips. Amber graduations are acid and alkali resistant.  FORTUNA® burettes are calibrated to deliver ('EX'). Please remember the waiting time of 30 sec.

Stopcocks:  FORTUNA® burettes are offered with straight stopcocks and glass or PTFE keys.  FORTUNA® Automatic Burettes, Pellet pattern are offered with lateral stopcock and glass or PTFE keys, with or without intermediate stopcock.

Precision Burette Tips: During the production of burettes we take great care with the burette tips so that the user can always titrate with the same size of drop.

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
volac FORTUNA® Burette, with straight stopcock and glass key

Class AS

conformity approved,
borosilicate glass 3.3

1 00802G, 04G

Schellbach,
blue or amber graduation

RB9803/NC

Schellbach, blue graduation

RB9803/WAC

Schellbach, blue graduation,
with single certificate

1 00806G

glass amber stain diffused,
white graduation

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volac FORTUNA®
Burette, with straight stopcock and PTFE key

Class AS,
conformity approved,
borosilicate glass 3.3

1 008..02P, 04P

Schellbach,
blue or amber graduation

RB9806/NC

Schellbach, blue graduation

RB9806/WAC

Schellbach, blue graduation,
with single certificate

1 008..06P

glass amber stain diffused,
white graduation

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volac FORTUNA®
Burette, Class USP, AS

with interchangeable
glass key in glass barrel

US9803/WAC

with certificate

with interchangeable
PTFE key in glass barrel

US9806/WAC

with certificate

volac FORTUNA®
Bürette, Klasse USP, AS

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FORTUNA® Burette, Bang pattern

Class AS,

conformity approved (10 ml),
borosilicate glass 3.3

1 030

with straight glass stopcock and glass
key in upper stopcock, Schellbach,
blue graduation

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⁽¹⁾ class AS, without conformity approval

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Packing unit: 1

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FORTUNA® Automatic Burette, Pellet pattern

Class AS,
without conformity approval,
automatic zero adjustment,
borosilicate glass 3.3

1 055..02

without intermediate stopcock, with
lateral PTFE stopcock, Schellbach,
blue graduation

1 055..06

without intermediate stopcock, with
lateral PTFE stopcock, amber stain
diffused, white graduation

Class AS,

conformity approved,
automatic zero adjustment,
borosilicate glass 3.3

1 057..02

with intermediate glass stopcock
and lateral glass stopcock, Schell-
bach, blue graduation

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
FORTUNA® Automatic Burette, Pellet pattern

Class AS,
without conformity approval,
automatic zero adjustment,
borosilicate glass 3.3

1 059..02

with intermediate PTFE stopcock and
lateral PTFE stopcock, Schellbach,
blue graduation

1 059..06

with intermediate PTFE stopcock and
lateral PTFE stopcock, amber stain
diffused, white graduation

Reservoir bottle for automatic burettes

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
Packing unit: single pack,
without reservoir

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FORTUNA® Automatic Burette, Dr. Schilling pattern

Class B

automatic zero
adjustment,
soda lime glass

1 068

Schellbach,
blue graduation

1 069

glass amber stain diffused,
white graduation

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
comes mounted on a polyethylene bottle, with a plastic foot

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Packing unit: 1

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FORTUNA® Graduated Pipettes

 FORTUNA® graduated pipettes are produced to the highest quality standards and are produced in Class AS. **(EN ISO 835).**

Reading:  FORTUNA® graduated pipettes are produced with ring marks. Reading of the volume at the lowest point of the meniscus.

Technical Details:

- graduated pipettes are calibrated to deliver
- Colour-Code-System (ISO1769)
- short practical delivery time, openings uniformly calibrated
- tips and suction tubes reinforced
- fire polished facettes, bevels surface tempered
- amber graduation acid and alkali resistant
- indication of production batch (for conformity approved pipettes)
- soda lime glass

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volac FORTUNA® Graduated Pipette

total delivery, **Class AS**,
conformity approved,
clear glass, with ring marks

1 1144..02,..04

blue or amber graduation

R400 - Type 3

blue graduated, zero at top,
capacity at outlet, total
delivery, non blowout

R385/WAC - Type 3

blue graduated, zero at top,
capacity at outlet, total delivery,
non blowout, with single certificate

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Original FORTUNA® Graduated Pipette with suction piston

clear glass, **Class A**, clear calibrated
cylinders, without certificate

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* additional to standard, without conformity approval

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volac FORTUNA® Graduated Pipette

USP Class AS, blue graduation. Each pipette is individually numbered and fully traceable satisfying USP/ISO/BS/DIN accuracies.

US385/WAC - Type 3

zero at top, capacity at outlet, total delivery, with single certificate

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Packing unit: 5

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FORTUNA® **Pipette filler**

the most common pipetting aid for all types of pipettes

The filler has three glass ball valves which can be controlled by finger pressure.

Functions: releasing the air, filling the pipette, delivery of the liquid

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Packing unit: 1

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FORTUNA® Bulb Pipettes

 FORTUNA® bulb pipettes are produced to the highest quality standards. They are offered in **class AS conformity approved** ('A' stands for the highest quality specifications, 'S' for swift delivery).

 FORTUNA® bulb pipettes are calibrated to deliver ('EX') and are supplied with 1 or 2 ring marks. **(EN ISO 648)**.

Reading: Reading of the volume at the lowest point of the meniscus.

Technical Details:

- Colour-Code-System (ISO 1769)
- short practical delivery time, openings uniformly calibrated
- tips and suction tubes reinforced
- fire polished facettes, bevels surface tempered
- amber graduation acid and alkali resistant
- indication of production batch
- soda lime glass

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volac FORTUNA® Bulb Pipette

Class AS, with 1 mark,
clear glass, conformity approved,
with batch certificate

R371/NC
blue graduation

R371/WAC
blue graduation,
with single certificate

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volac FORTUNA® Bulb Pipette

Class AS, with 2 marks,
clear glass, conformity approved,
with batch certificate

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Packing unit: 5

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Original FORTUNA® Bulb Pipette with suction piston

clear glass, **Class A**,
clear calibrated
cylinders, with 1 mark,
without certificate

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volac FORTUNA® Bulb Pipette

USP Class AS, with 1 mark,
blue graduation

Each piece is individually
numbered and fully traceable
satisfying USP/ISO/DIN/BS
accuracies.

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volac FORTUNA® Cylinders

volac FORTUNA® graduated and mixing cylinders (borosilicate glass 3.3) are 'Made in Germany' under the highest quality standards. They are produced in **Class B** or **Class A** (conformity approved). (**EN ISO 4788**).

Reading: **volac** FORTUNA® graduated and mixing cylinders are calibrated to contain ('IN'), with short graduation marks (Class B) or ring marks (Class A). Reading of the volume is at the lowest point of the meniscus.

Technical Details:

- delivered with hexagonal base or with plastic feet (autoclavable up to 130 °C, Class B)
- with tolerance indications
- amber graduations acid and alkali resistant
- indication of production batch (for conformity controlled cylinders).

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volac FORTUNA® Graduated Cylinder

high form, **Class A**,
borosilicate glass 3.3,
hexagonal base

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Packing unit: 2 (*1)

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FORTUNA® Graduated Cylinder

low form, **Class B**,
borosilicate glass 3.3,
hexagonal base

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Packing unit: 2

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
volac FORTUNA® Mixing Cylinder

Class A, conformity approved,
borosilicate glass 3.3,
hexagonal base,
with PP-stopper

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volac FORTUNA® Cylinders

All **volac** FORTUNA® cylinders are made of borosilicate glass, have hexagonal feet and are printed in blue graduation.

Technical details:

- every spouted cylinder (graduated) is formed to ensure clean, non drip pouring
- every stoppered cylinder (mixing) has a tooled and ground neck fitted with a 'Polystop' PP interchangeable stopper
- all cylinders are individually calibrated and individually printed to ensure that the scale length exactly fits the calibration of each cylinder,
- the Class A cylinders are batch tested with a lot number, WAC cylinders are individually tested and numbered
- if a Works Certificate is required this will show the date of testing, test results at 3 points and the identification number of the cylinder.
- All Class A cylinders satisfy the requirements of International Standards (**EN ISO 4788**).

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volac FORTUNA® Graduated Cylinder

Acc. to ISO, Class A,
blue graduation,
borosilicate glass,
hexagonal base

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
volac FORTUNA® USP Graduated Cylinder

USP Class A,
blue graduation,
borosilicate glass,
hexagonal base

US263/WAC
with individual serial number
and single certificate

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
R264

volac FORTUNA® Mixing Cylinder

Class A,
conform to ISO,
blue graduation,
borosilicate glass,
hexagonal base,
with PP-stopper

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US264

volac USP FORTUNA® Mixing Cylinder

USP Class A,
blue graduation,
borosilicate glass,
hexagonal base,
with PP-stopper

US264/WAC

with individual serial number
and single certificate

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Packing unit: 2 (*1)

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FORTUNA® Volumetric Flasks

 FORTUNA® volumetric flasks are produced to the highest quality standards in borosilicate glass only. They are manufactured in class A, conformity approved. Volumetric flasks are available with rim, with ground neck (TS), with PP- or glass stopper. **(EN ISO 1042).**

Reading: FORTUNA® volumetric flasks are calibrated to contain ('IN') and are offered with coloured ring mark. Reading of the volume is at the lowest point of the meniscus.

Technical Details:

- borosilicate glass 3.3
- accuracy class A
- with tolerance indications
- amber ring marks acid and alkali resistant
- indication of production batch

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1 514

volac FORTUNA® Volumetric Flask

Class A, conformity approved,
borosilicate glass 3.3

1 500

amber graduation,
with rim

1 514..02F

blue graduation,
with TS, without stopper

1 514..04F

amber graduation,
with TS, without stopper

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volac FORTUNA® Volumetric Flask

Class A, conformity approved,
borosilicate glass 3.3

1 513..02F

blue graduation,
with glass stopper

1 513..04F

amber graduation,
with glass stopper

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1.515

volac FORTUNA® Volumetric Flask

Class A, amber stain diffused,
conformity approved,
borosilicate glass 3.3

1 515

white graduation,
with PP-stopper

1 516

white graduation,
with glass stopper

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Packing unit: 5 (*2)

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FORTUNA® Volumetric Flask, Kohlrausch pattern

Class A, without conformity approval, borosilicate glass

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Packing unit: 2

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volac FORTUNA® Volumetric Flasks

We offer a range of flasks from 1 ml to 5000 ml capacity, all manufactured to class A (ISO) standards, with permanent marking and tooled ground neck, fitted with interchangeable polypropylene stopper. All flasks are made of borosilicate glass and are available in clear or amber glass (**EN ISO 1042**).

Class A Accuracy

This range has been enlarged to cover all the former class B capacities. All flasks are batch tested for accuracy and are supplied with a batch test certificate and a LOT number.

Class A Accuracy, Type WAC

Each flask is individually tested and carries its own individual number, with works certificate of calibration.

Class A according to USP

Each flask is individually tested and carries its own individual number. A Works Certificate of calibration can be provided at extra cost.

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volac FORTUNA® Volumetric Flask

Class A, conform to ISO,
blue graduation, borosilicate
glass, with Polystop®

RB258

with batch certificate

RB258/WAC

individually calibrated,
serially numbered with
single certificate

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RB258/AM

volac FORTUNA® Volumetric Flask

acc. to ISO, Class A,
amber stain diffused,
white graduation,
borosilicate glass,
with Polystop®

RB258/AM

with batch certificate

RB258/AM/WAC

individually calibrated, serially
numbered, with single certificate

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Capacity Inhalt	TS NS	white/weiß	white/weiß
		Art.-No. Art.-Nr.	Art.-No. Art.-Nr.
5 ml	10/19	RB258/AM/A/5	RB258/AM/WAC/A/5
10 ml	10/19	RB258/AM/B/5	RB258/AM/WAC/B/5
15 ml	10/19	RB258/AM/C/5	RB258/AM/WAC/C/5
20 ml	10/19	RB258/AM/D/5	RB258/AM/WAC/D/5
25 ml	10/19	RB258/AM/E/5	RB258/AM/WAC/E/5
50 ml	12/21	RB258/AM/F/5	RB258/AM/WAC/F/5
100 ml	12/21	RB258/AM/G/5	RB258/AM/WAC/G/5
200 ml	14/23	RB258/AM/H/2*	RB258/AM/WAC/H/2*
250 ml	14/23	RB258/AM/I/2*	RB258/AM/WAC/I/2*
500 ml	19/26	RB258/AM/J/2*	RB258/AM/WAC/J/2*
1000 ml	24/29	RB258/AM/K/2*	RB258/AM/WAC/K/2*

Packing unit: 5 (*2)

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volac FORTUNA® USP Volumetric Flask

USP Class A, blue graduation, each flask is individually tested and serially numbered and fully traceable satisfying USP/ISO accuracies, with Polystop®

US258/WAC

with individual serial number and single certificate

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US258/AM

volac FORTUNA® Volumetric Flask

USP Class A, amber stain diffused, white graduation, each flask is individually tested, serially numbered and fully traceable satisfying USP/ISO accuracies, with Polystop®

US258/AM/WAC
with individual serial number and single certificate

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white/weiß

Capacity Inhalt	TS NS	Art.-No. Art.-Nr.
5 ml	10/19	US258/AM/WAC/A/5
10 ml	10/19	US258/AM/WAC/B/5
15 ml	10/19	US258/AM/WAC/C/5
20 ml	10/19	US258/AM/WAC/D/5
25 ml	10/19	US258/AM/WAC/E/5
50 ml	12/21	US258/AM/WAC/F/5
100 ml	12/21	US258/AM/WAC/G/5
200 ml	14/23	US258/AM/WAC/H/2*
250 ml	14/23	US258/AM/WAC/I/2*
500 ml	19/26	US258/AM/WAC/J/2*
1000 ml	24/29	US258/AM/WAC/K/2*

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volac FORTUNA® Trapezoidal Flask

Class A, blue graduation,
borosilicate glass, with
PP-stopper TS 7/16

RB257

clear glass, with batch certificate

RB257/WAC

clear glass, with single certificate

RB257/AM

amber stain diffused, with batch
certificate

RB257/AM/WAC

amber stain diffused, with single
certificate

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Packing unit: 5

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FORTUNA® Precision Gas Syringes (Catheter Syringes)

It is well known that gas volumetric analyses are indispensable to a better understanding of the fundamental principles of chemistry. Up to now, however, they have always been rather complicated as regards to the instruments involved (gas burette, Hempel's absorption pipette, explosion or slow burning gas pipette) and seen in the light of the difficulties facing the chemistry teacher, it was difficult to be sure of success with them.

The FORTUNA® precision gas syringe (catheter syringe) not only eliminates all these difficulties but makes tests so much quicker and easier to carry out so that it is already being used in both schools and universities in many countries. Its popularity is largely due to the simplicity of design and the exclusion of confining liquid, thus eliminating solution errors and enabling tests to be carried out with gases highly soluble in water.

Technical Details:

- clear neutral glass
- calibrated cylinders
- acid and alkali resistant amber stain diffused graduation
- FORTUNA® gas syringes are only for technical use

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2 012

2 000

The quantitative combustion of gases in the glow wire bulb is very much simplified too. Some of the tests which are important to teaching are as follows:

- analysis of inhaled and exhaled air (determination of breathing quotients)
- indirect water analysis
- combustion of elements for atomic weight analysis
- quantitative gas combustion (H_2 , CO , NH_3 , hydrocarbons)
- photosynthesis of hydrochloric acid

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Iodine number flask, Sendtner pattern

2 120

with hexagonal stopper,
borosilicate glass

2 120

2 152

FORTUNA® Graduated Cylinder for compressed volume

2 152

Class B, high shape, blue graduation,
borosilicate glass, round base,

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D812

volac Glass Pasteur Pipettes Made in the UK

All **volac** Glass Pasteur Pipettes are precision made from carefully selected tubing on fully automatic machines ensuring total uniformity of production. Tips and tapers are centred and can be checked by rolling the pipettes on a flat surface. Length and taper of the tips are carefully controlled. Constrictions for plugging are neatly formed at a uniform distance from the fire polished end. Constrictions are heat formed which avoids weaknesses or contamination of the glass found when moulding or tooling processes are used (ISO7712).

volac Glass Pasteur Pipettes are made with a slightly heavier wall thickness to provide the best quality product allied to greater strength to give lower breakages during transit and reduced risk of breakage in use or when fitting PVC teats. **volac** Glass Pasteur Pipettes are available either plain, plugged with non-absorbent roving or pre-plugged and sterilised in ethylene oxide (EtO).

Technical Details:

- available in standard sizes 150 mm, 230 mm, plain or plugged, sterilised or non-sterile
- made of soda lime glass, soda lime glass with cerium or borosilicate glass
- '50 dropper' pipettes also available
- superior wall thickness as standard

Technical Details

D810, D812 **Standard Pasteur Pipettes**

Top Body outer ϕ	6.5 mm - 8.0 mm
Wall thickness	0.5 mm - 0.6 mm
Jet outer ϕ	1.2 mm - 1.8 mm
Jet inner ϕ	0.7 mm - 1.5 mm

Pasteurpipetten aus Glas Made in the UK

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
plugged, non-sterile,
soda lime glass

D810/PP

plugged, non-sterile,
soda lime glass with cerium

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50-Dropper, unplugged or plugged, soda lime glass

D813

Teat for Pasteur Pipettes made of PVC, smooth for safety and durability

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04. Syringes · Spritzen

FORTUNA® OPTIMA® All Glass Syringe

Borosilicate glass, interchangeable,
only for technical purposes

7 102

clear glass, with glass tip,
Luer, amber graduation

7 140

clear glass, with Luer-Lock
tip, amber graduation

All Glass Syringe, model ‚Tuberculin‘

Borosilicate glass, long type,
amber graduation, only for
technical purposes

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7 200

SANITEX Glass and Metal Syringe

Borosilicate glass, with steel piston, interchangeable, amber graduation, for technical use only

7 200

with Record tip

7 202

with Luer tip

7 240

with Luer-Lock tip

SANITEX Glass and Metal Syringe, model ‚Tuberculin‘

Borosilicate glass, long type, with steel piston, interchangeable, amber graduation, for technical use only

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
04. Syringes · Spritzen

7 150

FORTUNA® Glass Syringes for DGA oil-gas-analysis (without accessories)

Borosilicate glass,
glass piston,
amber graduation,
for technical use only

7 150

with Luer-Lock tip

Accessories for DGA Syringes

7 15101

Un-vented Male Plug Cap for
DGA-syringe

7 15102

3-way-Stopcock for DGA-syringe,
made of Nylon (price per m)

7 15103

Tube für DGA-syringe, made of
Tygon®, with Connector

7 15141 up to -51

Padded Carrying Case for DGA-
syringe

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04. Syringes · Spritzen

**Single Use Syringe
HENKE-JECT® 2-part,
without needle, steril,
CE**

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04. Syringes · Spritzen

7 410

**Single Use Syringe
HENKE-JECT®**, 2-part,
with needle, sterile, CE

7 470

with Luer tip

**Single Use Syringe, HENKE-
JECT®**, 2-part, without
needle, not sterile, CE,
bulk pack

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Packing unit: 100

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04. Syringes · Spritzen

Single Use Syringes
HENKE-JECT®,
3-part, without needle,
sterile, CE

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7 451

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7 460

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Single Use Syringes
HENKE-JECT®,
3-part, with by-
packed needle,
sterile, CE

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7 455

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Single Use Needles,
HENKE-JECT®,
ISO 7864

standard, sterile,
non-corroding steel, CE

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Packing unit: 100

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Poulten & Graf Laboratory Products for your needs

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05. Laboratory Products · Laborprodukte

For all the following products we guarantee the quality of raw materials and precise performance of the finished products. The user of glass and plastic labware can be confident that they can count on laboratory equipment from Poulten & Graf.

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
Borosilicate glass

8 080

with interchangeable plastic knob
lid and porcelain plate

8 100

with tubes (outlet) and stopcock TS
24/29 in lid and porcelain plate

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6 946

Glass stopper, DIN 12252

with hexagonal plate,
hollow, borosilicate glass

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
9 700

101305

Laboratory Bottle with screw thread, EN ISO 4796-1

8 780

clear, white graduation, autoclavable, borosilicate glass, with cap

8 781

amber, white graduation, autoclavable, borosilicate glass, with cap

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Amber Glass Bottle with screw thread

9 700

non coated, 250 - 1000 ml: square

9 710

plastic coated, square
1000 ml: also available round shape

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Storage Bottle, PTFE

101305

with screw cap (suitable for FORTUNA UNIVERSAL® HF dispenser)

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Technische Informationen

What is understood by quality?

In general usage the term 'Quality' has a number of meanings. Glass products of first class quality might be called this because the finish is good, the raw material used is good, the packaging is good, the appearance is good or the accuracy is good. Quality is defined as having all those characteristics and features of a product, which refer to its suitability for the fulfilment of given requirements, that is 'fit for purpose'. So Quality is related to the need and the specification set.

The Quality Management of Poulten & Graf

Product quality is directly dependent on the quality of process. Quality assurance must begin in the very first phase of the product creation process and cover all following phases. This global view of quality assurance, involving all aspects of production, design, testing, packing and conformity with both internal standards and externally set specifications, is the basis of our quality management. This is reflected in the approvals of ISO 9001/2015, which we have achieved and which we maintain.

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06. Technical Support · Technische Details

Certification or attestation of conformance

Since January 01, 2015, new German Metrology and Verification Regulations have come into force. They replace the previous German Weights and Measures Regulations. For Class A and Class AS Volumetric Glassware which are produced in Germany, a new conformity sign is being used. The H-sign is now replaced by the **DE-M**-sign.

As before, the manufacturer is allowed to certify the conformity on basis of the PTB-regulations (National Metrology Institute). Government institutions will ensure that the regulations are being followed. Therefore, nothing has changed in the production and testing methods. If required, a manufacturer's declaration of conformity can be enclosed with class A and AS volumetric products.

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Works Certification of Accuracy of Volumetric Glassware

All  FORTUNA® Class A/AS glassware is conformity approved. Test records are kept in house with the manufacturing records. However individual Works Certificates of Accuracy can be provided for each piece if required; each piece of glass has its own individual number engraved on it and the Certificate supplied will have the corresponding number, the date and the relevant test results.

Such Works Certificates of Accuracy (look for the designation WAC in the product codes) should not be confused with Certificates of Conformity which refer to batch testing only.

Poulten & Graf Ltd. 1 Alfreds Way, IG11 OAS Barking, UK	
Qualitätszertifikat Certificate of Performance Certificat de Qualité - Certificado de Calidad	
Seriennummer: Series - Série - Serie:	V000153
Artikelnummer: Code - Número del Artículo:	US258/WAC/G/5
Artikelbezeichnung: Article - Artículo:	VOLAC Flasks, 100ml, TS12, USP, ind. Cert.
Nennvolumen: Nominal volume - Volume nominal - Volumen nominal:	100 ml
Toleranz in % und ml: Tolerance in % and ml / Tolérance de % et ml	0,080% 0,080 ml
Norm: Standard - Norma:	EN ISO 1042 / ASTM 288
Hiermit bestätigen wir, dass die oben aufgeführten Messgeräte nach den Vorgaben der Eichordnung einzeln geprüft sind. Alle ermittelten Messwerte liegen innerhalb der angegebenen Toleranzgrenzen.	
Gemessene Auslaufzeit: measured time of flow - écoulement mesure	0
Datum: Date - Fecha:	10.03.2016 11:28:15
Prüfer: Operator - Vérificateur - Comprobador:	Passed Weis

The system for LOT numbers on FORTUNA® volumetric glassware

All  FORTUNA® conformity approved volumetric glassware is delivered with the FORTUNA® Batch Certificate of Performance. A regular control of all testing equipment is required by DIN EN ISO 9001. The requirement for traceability is met by the batch certificate and with the easily identifiable LOT numbers (e.g. LOT B00008) on  FORTUNA® volumetric glassware, which is produced in Wertheim.

Beside the conformity sign all  FORTUNA® volumetric glassware is marked with an individual LOT number that allows an exact classification back to the batch certificate.

If your batch certificate is lost, you may ask for a copy by informing us of the product details and the LOT number or request it online.

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A00184

To define the different classes (Class A, Class AS, Class B, etc.) for glass or 'performance' for pipetors and dispensers the following terms are used: 'accuracy' and 'coefficient of variation' (precision). Sometimes we also speak about 'tolerances' which are specified in the ISO-standards. This means the permitted deviation of the measured values from the nominal volume. By accuracy we mean corresponding to the average value of a series of measurements adjusted for temperature and pressure with the desired value. Precision is a measure of the variation from the nominal value achieved through a series of measurements at one setting. The coefficient of variation is the % value of the precision.

Measurement of volumes is part of the routine work in a laboratory. The most common volumetric glassware such as bulb- and graduated pipettes, burettes, volumetric flasks and graduated cylinders form part of the basic equipment in a laboratory. Inscriptions and graduations on volumetric glassware should be easy to read and acid- and alkali-resistant. The complete range of  FORTUNA volumetric glassware fulfils these requirements.

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
Calibration of Volumetric Glassware

- **EX** calibrated to deliver. The delivered liquid corresponds with the volume indication. The liquid that remains on the wall or in the tip of the glassware is taken into consideration.
- **EX + s** calibrated to deliver after a waiting time (for pipettes 5 seconds, for burettes 30 seconds). Keeping to the waiting time is important to avoid mistakes in measurement.
- **IN** calibrated to contain. The contained liquid corresponds with the volume indication. Example: flasks and cylinders.
- **20°C** The calibration has been done at 20°C. Volumetric glassware changes volume about 0,1 % each degree Celsius.

Precision Classes

- **A** tolerance according to ISO standard. Class A is the most precise class. Glassware of class A is suitable for official calibration, this means it can be tested by either the German or the UK Weights and Measures Offices. This is useful for control of measuring instruments according to DIN ISO 9001.
- **AS** as class A, but for pipettes and burettes with fast delivery. The waiting time is much shorter than for class A.
- **B** the tolerances of class B are about twice those of class A and AS. Volumetric glassware of class B is suitable for routine work in the laboratory.

Meniscus

Value is read at the lowest point of the liquid surface, which means when the lowest point of the meniscus touches the upper line of the graduation mark (see picture 1).

For burettes with Schellbach stripes the reading has to be made at the level where the two wedge-shaped points are touching (see picture 2). When reading the scales it is important that one's eyes are at height of the fluid level. Otherwise it results in reading errors (parallax errors).

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
A Guide to calibrate Micro Pipettors and Bottle Top Dispensers (compiled from BS ISO 8655: Part 1-6)

Test Lab Requirement

- Lab should ideally be air conditioned.
- Temperature and air pressure should be stable.
- Air movement should be minimised.
- Weighing/balance surface should be free of vibration.

Test Lab Equipment

- Air pressure barometer
- Certified Class A thermometer
- Accurate balance (5 decimal places)
properly calibrated to National Standards
- Distilled water at stable room temperature

Pipetting Practice for Accurate Results

Use high quality tips preferably of a brand to match the pipettor. Ensure tips do not leak (hold pipettor vertically, tip filled with liquid, without dripping for at least 20 secs). Rinse and pre-wet each tip with distilled water at least 3 times before use. Fill and dispense slowly and smoothly. Fill with pipettor tip not more than 3-4 mm below liquid surface and pipettor held vertical (fig. 1) Dispense with pipettor tip against wall of container and pipettor held at 45° to the vertical and draw tip up the container wall one or two seconds after second blow out (fig. 2).

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

2

Pipetting - Balance Readings

Tare a container on the balance. Using a wetted tip, pipette/dispense into the tared container on the balance. Read weight of dispensed water. Read temperature of distilled water. Read the air pressure. For setting calibration test 4 times at each volume, for checking/certifying calibration test minimum 10 times and for testing coefficient of variation test 30 times.

Calculation of Accuracy

- Calculate mean weight \bar{W} from the number n of individual weighings W_i summed to give $\sum W_i$ using the equation:

$$\bar{W} = \frac{\sum W_i}{n}$$

- Calculate mean volume V_t of liquid samples at test temperature t :

$$\bar{V}_t = \bar{W} \cdot Z$$

where Z is the conversion factor for the density of water at given temperatures and air pressure - see table (page 252)

N.B. If an electronic balance is used, the effects of evaporation are minimal and can be ignored; if a mechanical balance is used mean evaporation loss \bar{e} needs to be estimated and added to \bar{W} .

- Calculate inaccuracy in percent of the mean \bar{A}_t at test temperature t using V_t mean volume and V_o nominal volume by equation:

$$A_t = \bar{V}_t - V_o \cdot \frac{100}{V_o}$$

Index of Symbols

\bar{W} mean weight

W_i results of individual weighings

n number of weighings

t water temp. of test

\bar{V}_t mean volume at test temperature

V_o nominal volume

\bar{e} mean evaporation loss

Z Conversion factor (see table page 254)

\bar{A}_t Accuracy (inacc.)

as percent of mean

S Standard deviation of weighings (or volumes)

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06. Technical Support · Technische Details

Calculation of Precision (Repeatability) - Coefficient of Variation

- Calculate standard deviation **S** by the equation:

$$S = \sqrt{\frac{\sum (W_i - \bar{W})^2}{n-1}}$$

(n should be not less than 10 and should ideally be 30)

- Calculate coefficient of variation **CV**:

$$CV = \frac{S}{\bar{W}} \cdot 100$$

If the balance is attached to computing equipment \bar{W}_i and **W** may be replaced by equivalent volumes.

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Calibration of Volumetric Glassware (see BS 6696)

Test Lab Requirement

- Full temperature control (air conditioning)
- Atmospheric Pressure Barometer
- Certified Thermometer
- Accurate calibrated weighing balance

Cleaning and Drying

Wash, clean and dry glassware thoroughly. Ensure that a satisfactory method of cleaning is employed, since any contamination will effect the meniscus shape/level.

During calibration it is critical that all internal surfaces are clean and dry to avoid moisture droplets or contaminants.

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

Timing Delivery

Delivery glassware should be timed before calibration, to ensure that total delivery time is within tolerances specified by the relevant BS and ISO standards.

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Weigh and Tare


With delivery glassware such as pipettes and burettes, weigh and tare a suitable sized weighing vessel.

Glassware should be maintained at a controlled room temperature (air conditioning).

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
Fill with purified water a few millimeters above the graduation mark. Pipettes and burettes should be filled from below, i.e. through the tip. If burettes are filled from the top for added convenience, take care to avoid wetting the walls above the upper graduation mark and check that no air is trapped in the stopcock. For pipettes and burettes control the flow of excess liquid from the tip to allow the meniscus to settle exactly on the graduation mark. For flasks and cylinders remove excess water to match meniscus and graduation mark exactly using a pipette and bulb.

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
When setting the meniscus ensure that the glassware is vertical. Ensure parallax error is avoided by viewing the meniscus from an eyelevel height against a white background.

In addition a black shade should be secured either around the article to be tested or placed behind it. The shade should be below the setting height (approx 1 mm lower). Ensure the meniscus edge is smooth, straight and level.

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06. Technical Support · Technische Details

For graduated/bulb pipettes and burettes, dispense water into a tared weighing vessel and weigh. For flasks and cylinders weigh the tared and filled glassware. Record the result from the balance. Use a certified Class A thermometer to measure the water temperature. Add the water temperature correction figure to the net weight obtained.

Measure the barometer reading and air temperature and add (or subtract if minus) the correction figure.

This produces the final gravimetric result which has already been converted to the volumetric figure in millilitres.

All of the procedure described must be followed to ensure that consistently accurate results are obtained.

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Environment and safety

Current discussions on environmental protection, dangerous materials etc. have made laboratory technicians more sensitive in the selection of laboratory instruments. Safety aspects remain the primary concern while handling dangerous or hazardous materials. Laboratory equipment must be extremely safe and reliable. All regulations must be observed to ensure protection of both humans and the environment.

We offer a repair service for a number of products:

- bottle top dispenser POLYFIX®, OPTIFIX®,
 - electronic burette OPTIFIX® titrier,
 - automatic dosing system OPTIMAT® 2 and 3
- Please remember that when returning any instrument for repair it must be accompanied by a written certificate of decontamination.

When disposal of waste is involved we ask you to consider the necessary regulations. These apply also to disposal of single use products.

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06. Technical Support · Technische Details

Important safety instructions:

1. Laboratory instruments and equipment should be checked on purchase as being in good condition and working order before use.
2. Laboratory instruments should be checked for damage as part of a daily routine to avoid risk of injuries such as cuts, burns and infections.
3. Laboratory glassware which shows any sign of damage should be disposed of in a safe manner.
4. Sudden changes in temperature should be avoided while handling glassware.
5. Glass equipment should be inspected thoroughly before using in a vacuum or under pressure.
6. Sudden changes in pressure should be avoided. Vessels should be vented smoothly when used in a vacuum. Flat bottomed glassware should not be used under pressure or in a vacuum.

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06. Technical Support · Technische Details

Cleaning of Laboratory Glassware

Laboratory glassware should be washed before first use.

Cleaning by hand

The common cleaning method is to wipe and rub the glassware with a cloth or sponge soaked in a cleaning solution. Abrasive cleaners and sponges should not be used because they can damage the surface of the glass.

Machine cleaning

Cleaning in a machine is more gentle for glassware than soaking it. Glass comes into contact with the dishwashing liquid for relatively short periods when it is being sprayed onto the surface of the glass.

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Cleaning of Bottle Top Dispenser

To ensure proper operation and life of the OPTIFIX® dispensers, we recommend the following cleaning procedure to provide you with years of trouble free service.

- Prepare a standard cleaning solution and hot water.
- Remove the piston from the cylinder via the volume setting knob. Remove discharge tube, intake tube and adapter.
- Put all parts into the cleaning solution for 2-4 hours or overnight to allow a good soaking.
- After soaking, use a clean glassware brush to remove any further dirt if necessary.
- Rinse all dispenser parts with distilled water and air dry.
- Reassemble the unit and put the OPTIFIX dispenser back into service.

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06. Technical Support · Technische Details

Disinfection

Laboratory equipment that has been in contact with infectious material should be disinfected before use, for the protection of laboratory workers. It can be washed by hand with a cleaner-disinfectant. Further it can be cleaned by using physical or thermal methods (for 10 minutes at 93°C as required by the BGA). If necessary the glassware can be steam sterilized.

Steam sterilization

Steam sterilization is 'the destruction or irreversible inactivation of all microorganisms of 120°C and 2 bar'. At a sterilization temperature of 121°C the glassware should be sterilized for 20 minutes minimum effective application time. Laboratory equipment should be always cleaned carefully before steam sterilization, otherwise any soiling will bake on during the steam sterilization. If there are any chemicals on the surface of the glassware during the sterilization they can damage the glass due to the high temperatures. Containers should be opened during the sterilization to avoid a build-up of pressure. The steam has to have unrestricted access to all contaminated points to ensure effective steam sterilization.

Raw materials we use

Soda Lime glass finds its use in the foodstuffs industry as bottle and preserve glass as well as in the laboratory industry. The chemical composition of this glass consists of 71-75 % sand (SiO_2), 12-16 % soda and 10-15 % lime. Soda lime glass is subject to acute stress under heat shock and has a relatively high coefficient of expansion. It is particularly resistant to alkali liquids.

Borosilicate glass differs from soda lime glass by using a higher proportion of sand (SiO_2). The composition of this glass contains 70-80% silicon dioxide (SiO_2), 7-13% boron trioxide (B_2O_3), 4-8% Na_2O and K_2O , as well as 2-7% alumina oxide (Al_2O_3). Glass with such composition is characterised by a high resistance against acidic attack, low thermal expansion and excellent temperature stability. This glass is particularly suitable for laboratories and chemical use where contact with aggressive materials is likely to be prolonged.

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06. Technical Support · Technische Details

The POLYFIX® Dispenser is made of Polypropylene and glass.

The OPTIFIX® dispenser is made of PTFE and glass, while the SMART pipettor is made of Polypropylene and Polycarbonate.

PP - Polypropylene

When polymerizing, propylene develops a product which contains methyl side groups of monomers of the paraffin chain. The volume requirement due to the regular methyl side chains is relatively large. For this reason polypropylene is the plastic with the lowest specific weight in compact form. It is similar to polyethylene, and may be used instead of polyethylene, where it particularly depends on good heat resistance, high impact strength, inherent stability during stress or particularly low specific weight. Polypropylene also has a good resistance to many inorganic reagents as well as the less aggressive organic solvents.

PTFE - Polytetrafluorethylene

PTFE has no thermoplastic characteristics although it consists of linear, chemical chain molecules which are not linked. It is not possible to handle or form PTFE in the usual thermoplastic way. Nevertheless PTFE is a thermoplast because it can be sintered. For the production of shaped parts from PTFE the polymer is cold pressed as powder under high pressure into the desired form and sintered at temperatures around 380°C. PTFE is attacked only by elemental fluorine and chlorine tri fluoride at higher temperatures and pressures as well as by melting alkali metals. Otherwise it is resistant to all chemicals. A further advantage of polytetrafluorethylene: it has the lowest coefficient of friction of all solids. However above 350°C PTFE burns in air to form noxious gases.

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Poulten & Graf GmbH – Conditions of Sale

1. General conditions

Our General Terms of Business are binding for all deliveries and services. We will not accept general terms of business of our contract partners, which differ from our terms of sale, even without explicit objection. Additional agreements have to be confirmed by us in writing. Only German Law and German Text will be applicable. Place of performance is Wertheim. The place for jurisdiction, also for bill transactions, is Mosbach. All important data necessary for business transactions may be stored on computer. All information, which is not public knowledge and which results from the business relationship with our company, must be kept secret from third parties.

2. Offers

Our offers are always without engagement. In order to be valid, all transactions, orders and other agreements require our confirmation in writing. Our written order confirmation or agreement shall determine the nature and scope of the delivery. We reserve the right to make technical alterations to the offered products.

3. Risk, conditions of delivery and non-acceptance of delivery

The risk of ownership will be transferred to the buyer, as soon as the goods leave our company even in the event of delivery carriage-paid. Complaints can only be taken into account if submitted within 8 days of receipt of the goods. Delivered goods will only be accepted for return in their original packing. Illustrations and brochures etc. are not binding. Measurements and utilisation data are only to be regarded as approximate and are not binding.

Deliveries will be made ex factory, packing and insurance excluded, if no other agreement exists. Packing material of any kind may be returned to us within the range of legal regulations applicable to us. We will not pay for the cost of return of packing materials. In the event that dispatch or delivery is delayed at the request of the customer or for reasons for which he is responsible, the risk shall pass to the customer for the period of the delay, however at the request of the customer, we will arrange insurance cover that he requires at his expense.

Partial deliveries are permissible. Surplus quantities within the normal range (up to 15 %) are permissible and must be accepted and paid for by the customer. In case of delayed acceptance of the goods, we may sell the items at the customer's expense to a third party or we may ask for compensation.

4. Delivery times and Unforeseen Circumstances

Indicated delivery times are without obligation, unless a binding delivery time has been fixed in writing. Deliveries may be delayed due to Act of God or strikes, interruption of work, shortage of raw materials or failure to deliver by our suppliers. Liability for delayed delivery will be limited to proven, gross negligence, caused by us.

5. Special productions

Products, which are manufactured according to drawings, samples or instructions of the customer, cannot be returned. This includes in particular individual software for computers, delivered by us. All kinds of software delivered are subject to copyright and may be neither copied, modified nor presented in public without our prior written permission. Additional quantities within the normal range (up to 15 %) are permissible and must be

accepted and paid for by the customer, unless exact quantities have been agreed upon in writing. The customer is responsible for ensuring that products manufactured according to his instructions, do not interfere with any trade right, copyright or other rights of a third party. The customer must compensate all damages resulting from any violation of such rights.

6. Payment conditions, prices, delayed payment

Products will be invoiced in EURO at the prices in effect on the day of delivery. Prices are ex factory excluding packaging and transportation insurance with immediate payment, unless otherwise provided. We reserve the right to include additional payment conditions for new customers. The customer may only set off payment against counter-claims, which are undisputed or legally determined. In case of delayed payment, we reserve the right to charge the customer with all costs as well as interest at 3 % above the respective EURIBOR 3 Month interest rate.

7. Warranty Claims

Warranty claims shall be recognised, if the defect has been reported to us in writing within 12 months of the delivery date. To the exclusion of further claims – subject to section 8 – we will provide a warranty as follows:

7.1 We shall have the option to either replace or to carry out repairs free of charge on all those parts which prove to be defective due to circumstances which occurred prior to the transfer of risk. The discovery of such defects must be reported to us immediately in writing. Parts replaced will become our property.

7.2 If we have to supply according to drawings, specifications, samples etc. provided by the customer, then the customer carries the risk of suitability for intended purpose. The deciding factor for the condition of the goods conforming to the agreement is the point in time of the transfer of risk according to section 3.

7.3 If a final inspection of the goods has been agreed or an inspection of the first sample, then customer's complaints will not be recognised later concerning defects which the customer should have noticed at the final inspection or inspection of the first sample if it had been carried out with due care and attention.

7.4 We must be given an opportunity to confirm the defect the complaint is being made about. Rejected goods must be returned to us immediately on demand. We will bear the transport costs if the customer complaint is justified. If the customer does not comply with these obligations or makes changes to the rejected goods without our agreement, then he forfeits any claims under warranty of quality.

7.5 We will bear the costs for replacement of goods including transport and reasonable costs for disassembly and assembly incurred due to the rectification of defects or the replacement delivery, providing the complaint proves to be justified.

7.6 Within the framework of legal provisions, the customer has the right to withdraw from the contract if, taking the legal exceptions into consideration, we do not successfully make use of a reasonable period of time given to us for rectification of defects or a replacement delivery due to a defect. If the defect is not significant, then the customer only has the right to claim a reduction of the contractual price. In all other cases the right to a reduction of the agreed price is excluded.

7.7 No warranty will be given in the following cases in particular: Unsuitable or inappropriate use, faulty installation or start-up by the customer or a third party, normal wear

07. Terms & Conditions · AGBs

and tear, incorrect or negligent treatment, maintenance not carried out correctly, chemical, electrochemical or electrical influences providing we are not liable for them.

7.8 If the customer or a third party reworks the goods incorrectly, then we are not liable for the consequences arising from this. The same applies to changes made to the delivered item without our prior agreement. 7.9 We only provide a warranty for accessories added on within the framework of the warranty provided by our supplier.

8. Liability

8.1 If we are responsible for the customer not being able to use the delivery item as stipulated in the contract as a consequence of suggestions made and consultations carried out before or after the contract was concluded which were omitted or carried out incorrectly or through infringing upon other subsidiary obligations - in particular instructions for operating and maintaining the delivery item - then to the exclusion of further claims made by the customer, the stipulations in sections 7 and 8.2 apply correspondingly.

8.2 For damages, which have not been incurred by the delivered item itself, we only accept liability - whatever legal reasons are given - in the case of

- intent,
- gross negligence by executive employees,
- culpable injury of life, body, health,
- malicious silence with regard to a defect or where a defect's non-existence had been guaranteed
- defects of the delivery item for which liability must be accepted according to German Law on Product Liability for damages to persons or property with regards to privately used objects. In the case of culpable infringement of essential contractual obligations, we also accept liability for gross negligence of non-executive employees and in the case of slight negligence, in the latter case limited to reasonable, foreseeable damages typical to a contract.

Any further claims shall be excluded.

9. Reservation of proprietary rights

The delivered goods remain our property until complete payment. If the goods are resold or processed by the buyer, we are to be considered as manufacturer according to § 950 BGB and achieve the property of all intermediate or final products. The processing party is only the depositary.

For other items, which are not our property, we achieve the co-property of this new product in proportion to the goods reserved to the other components. The goods may only be sold within ordinary and regular business transactions and only, if claims arising from resale have not been assigned to a third party. The claims of the buyer resulting from resale are to be considered as assigned to us upon conclusion of the contract of sale. This is also valid, if our goods are combined or processed with other items. In this case, the assigned claims serve as our security only to the extent of the value of the merchandise sold under reservation of proprietary rights. The buyer is obliged to provide us with the names of third party debtors upon request and to inform them about the assignment. Amounts, which have been collected by him, have to be paid over to us, if our claims are due for payment. Distraints and hypothecations of the goods under reservation of proprietary rights or of the assigned claims are not permissible. The buyer has to inform us immediately about any access of a third party to the goods under reserve

of proprietary rights or to the claims assigned to us. In the event of sales to foreign countries, should the retention of title not be admitted with the same effect as under German Law, the goods involved will remain our property until payment of all our claims resulting from the contractual relationship brought about by sale of said goods. In case, however, this retention of title is not admitted with the same effect as under German Law, but if it is allowed to reserve other rights to said goods, we will be entitled to exercise all of these rights. The buyer will be obliged to co-operate in all measures that we take for the protection of our right of ownership or, in its place, of any other title to the goods.

10. Tools

Tools, moulds or other devices incorporating our know-how, which have been produced by us or handed over to the customer, remain our property, even if the customer has paid the costs for these in part or in full.

11. Statutory Limitation

All claims made by the customer - for whatever legal reasons - are subject to a limitation period of 12 months. For wilful or fraudulent behaviour and for claims based on the German Product Liability Law, the legal limitations apply.

12. Validity of the Agreement

Should one or more of the provisions of this Agreement be or become invalid for any reason, the provisions should be interpreted to ensure that the commercial aim of the original, invalid provision is upheld. The validity of the remaining provisions will remain unaffected.

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