

04. Syringes · Spritzen

7.200



SANITEX Glass and Metal Syringe

Borosilicate glass, with steel piston, interchangeable, amber graduation, for technical use only

7.200

with Record tip

7.202

with Luer tip

7.240

with Luer-Lock tip

SANITEX Glass and Metal Syringe, model ‚Tuberculin‘

Borosilicate glass, long type, with steel piston, interchangeable, amber graduation, for technical use only

7.260

with Record tip

7.262

with Luer tip

7.264

with Luer-Lock tip



SANITEX Recordspritze

Borosilikatglas, mit Stahlkolben, austauschbar, braun graduiert, nur für technische Anwendungen

7.200

mit Record-Konus

7.202

mit Luer-Konus

7.240

mit Luer-Lock-Konus

SANITEX Recordspritze, Modell ‚Tuberkulin‘

Borosilikatglas, lang, mit Stahlkolben, austauschbar, braun graduiert, nur für technische Anwendungen

7.260

mit Record-Konus

7.262

mit Luer-Konus

7.264

mit Luer-Lock Konus

		7.200	7.202	7.240
Capacity	grad./ml	Art.-No.	Art.-No.	Art.-No.
Inhalt	grad./ml	Art.-Nr.	Art.-No.	Art.-No.
2 ml	: 0.1	7.200-27	7.202-27	7.240-27
5 ml	: 0.2	7.200-33	7.202-33	7.240-33
10 ml	: 0.5	7.200-37	7.202-37	7.240-37
20 ml	: 1.0	7.200-41	7.202-41	7.240-41
30 ml	: 1.0	7.200-44	7.202-44	7.240-44

Packing unit: 1

Verpackungseinheit: 1

		7.260	7.262	7.264
Capacity	grad./ml	Art.-No.	Art.-No.	Art.-No.
Inhalt	grad./ml	Art.-Nr.	Art.-No.	Art.-No.
0,5 ml	: 0.01	7.260-17	7.262-17	7.264-17
1 ml	: 0.01	7.260-19	7.262-19	7.264-19

Packing unit: 1

Verpackungseinheit: 1