

Superior Quality!

volac Glass Pasteur Pipettes



**Soda Lime
Glass**

**NEW:
CERIUM SODA
Glass**

**Borosilicate
Glass**

All **volac** Glass Pasteur Pipettes are precision made from carefully selected tubing on fully automatic machines ensuring total uniformity of production. Tips and tapers are centred and can be checked by rolling the pipettes on a flat surface. Length and taper of the tips are carefully controlled. Constrictions for plugging are neatly formed at a uniform distance from the fire polished end. Constrictions are heat formed which avoids weaknesses or contamination of the glass found when moulding or tooling processes are used (ISO7712).

volac Glass Pasteur Pipettes are made with a slightly heavier wall thickness to provide the best quality product allied to greater strength to give lower breakages during transit and reduced risk of breakage in use or when fitting PVC teats.

Glass Pasteur Pipettes are available either plain, pre-plugged with non-absorbent roving or plugged and sterilised in Gamma irradiation.

Technical Details:

- available in standard sizes 150 mm and 230 mm, unplugged or plugged, sterilised or non-sterile
- made of soda lime glass, CERIUM SODA glass (NEW), borosilicate glass
- '50 dropper' pipettes also available
- not a medical device



Glass Pasteur Pipettes, disposable

Soda Lime Glass



	unplugged non-sterile soda lime glass	plugged non-sterile soda lime glass
Length	Art.-No.	Art.-No.
150 mm	D810	D810/PP
230 mm	D812	D812/PP

Packing unit: 1000 pcs (4 inner packs x 250 pcs)

NEW: CERIUM SODA Glass

	plugged sterile CERIUM SODA glass
Length	Art.-No.
150 mm	D810/VA/CERIUM
230 mm	D812/VA/CERIUM

Packing unit: 250 pcs. (10 x 25)

Borosilicate Glass

	unplugged non-sterile borosilicate glass	plugged non-sterile borosilicate glass
Length	Art.-No.	Art.-No.
150 mm	D810/BORO	D810/BORO/PP
230 mm	D812/BORO	D812/BORO/PP

Packing unit: 1000 pcs (4 inner packs x 250 pcs)

Soda Lime Glass

	50-Dropper unplugged soda lime glass	50-Dropper plugged soda lime glass
Length	Art.-No.	Art.-No.
150 mm	D897	D897/PP

Packing unit: 1000 (4 inner packs x 250)

Teat for Pasteur Pipettes

made of PVC, smooth for safety and durability

Art.-No.

D813

Packing unit: 100



Poulten & Graf GmbH

Karl-Carstens-Str. 10,
97877 Wertheim, Germany
Tel.: +49 9342 9229-0
www.poulten-graf.de
sales@poulten-graf.com

All data are subject to change - Copyright:
Poulten & Graf GmbH, GERMANY 2023/06

